

USER MANUAL_ELDERLY

Guardian Angel (EN)



CERTH
CENTRE FOR
RESEARCH & TECHNOLOGY
HELLAS



Introduction

Guardian Angel monitors vital signs and activity to detect critical situations. In case of a measurement out of range it can send notification alert with to the formal/informal caregiver.

How to charge the tablet

1


Connect the charging cable to the tablet and the power adaptor, as shown in the figure above. Place it in the wall socket
Charge every 2-3 days for 3 hours!!



How to charge the wristband

1

Detach the sensor from the bracelet



2

Insert the sensor into the cradle, with the 2 gold plated connectors inserted head-in



3

Insert the cradle into the power adapter of the tablet



Place it in the wall socket.
Charge every 2-3 days for 3 hours!!

For technical support please call 6972801970 & 2311257610 – Weekdays 11.00 – 14.00