

USER MANUAL_ELDERLY

Physical Activity Monitoring (EN)



CERTH

CENTRE FOR
RESEARCH & TECHNOLOGY
HELLAS



Introduction

Guardian Angel monitors vital signs and activity to detect critical situations. In case of a measurement out of range it can send notification alert with to the formal/informal care-giver.

How to charge the tablet

1


Connect the charging cable to the tablet and the power adaptor, as shown in the figure above. Place it in the wall socket
Charge every 2-3 days for 3 hours!!



How to charge the wristband

1

Detach the sensor from the bracelet



2

Insert the sensor into the cradle, with the 2 gold plated connectors inserted head-in



3

Insert the cradle into the power adapter of the tablet



Place it in the wall socket.
Charge every 2-3 days for 3 hours!!

For technical support please call 6972801970 & 2311257610 – Weekdays 11.00 – 14.00

Instructions for the User: How to use the Blood Pressure Monitor

1

Ensure that the Lock switch on the top of the monitor shows green colour



2

Wear the sleeve as shown below



3

Press the middle button with the label



4

The sleeve will inflate and you will see your blood pressure measurement and pulses on the screen.

5

Wait until the display shows OK and then press the middle button (of step 3) to turn the blood pressure monitor off.

How to replace the batteries

When the display shows the indication shown below, open the cover on the back of the monitor and replace the batteries.



Instructions for the User: How to use the Oximeter

1

Wear the oximeter on a finger and press the green button.



2

You will see on the screen your oxygen saturation and heart rate measurements

3

Wait until the oximeter screen turns off.

How to replace the batteries

When the display shows the indication shown on the left image, please open the lid and replace the battery.

